



FICHE DE DONNÉES DE PRODUIT

CA[®] Stripping Tools / CA[®] Strip Opener

Article		Matériau
REF 5342.1 – 5354.1	CA [®] Stripping Tools	A, B, C
REF 5326.1 + 5327.1	CA [®] Strip Opener	A, B, C

Matériau								
A	Acier inoxydable 1.4301							
	C %	Si %	Mn %	P %	S %	Cr %	Ni %	Autres %
	≤ 0,08	≤ 1,0	≤ 2,0	≤ 0,045	≤ 0,030	18,0 – 20,0	8,0 – 10,5	Se ≤ 0,15

B	Diamant
---	---------

C	Silicone
---	----------